

## CellSafe+ Biopsy Capsule

Please  
contact us for  
more information  
[cellpath.com](http://cellpath.com)

### Secure Biopsy Specimen Handling

The CellSafe+ Biopsy Capsule offers high specimen security. It is an easy-to-use alternative to the more traditional and time-consuming handling techniques such as tissue wrap and foam pads. The extra fine X-Ray transparent mesh also offers zero tissue damage and low reagent carry over.

The capsule is closed in one direction to fully encapsulate the tissue. During processing, CellSafe+ is contained in a standard cassette for processing in the normal way (for best results use with the System II Hex Cassette and lid).

The blue capsule aids contrast and visualisation of light coloured biopsies, normally experienced with breast tissue.

Part No	Description		Quantity
EBE-0201-02A	CellSafe+ Biopsy Capsule	White <input type="checkbox"/>	100
EBE-0201-03A	CellSafe+ Biopsy Capsule	White <input type="checkbox"/>	1000
EBE-0301-02A	CellSafe+ Biopsy Capsule	Blue <input checked="" type="checkbox"/>	100
EBE-0301-03A	CellSafe+ Biopsy Capsule	Blue <input checked="" type="checkbox"/>	1000

N.B. If using any cassette not manufactured by CellPath, compatibility tests must be carried out before routine work.

### V2 | Product UPGRADE

Firmer mesh is more robust  
providing added security



New catch feature gives a positive "click" sound and feel when closed, giving the user notification that the capsule is truly shut.

### KEY BENEFITS

- 1 Blue mesh aids contrast and visualisation of light coloured biopsies
- 2 X-Ray transparent for breast work
- 3 Extra fine mesh provides up to 20 times less carry over of solvents than foam pads
- 4 Zero tissue damage or artefacts associated with foam pads
- 5 Substantial time savings when compared to wrapping biopsy in tissue paper

00.784.285 V4