

The Diagnostic BioSystems **HighLighter**[™]

Automated Slide Stainer

The Diagnostic BioSystems HighLighter[™] is a Fully-Automated Slide Staining System

- Compact, benchtop, 36-slide system
- Designed for IHC application
- Fully automated from baking to counterstain
- Consistent and crisp staining
- Real time onboard DAB mixing
- Liquid sensors for accurate reagent handling
- Liquid sensors for waste disposal to prevent overflow
- Slide Data Matrix Code labeling and reading
- RFID reader for reagents vials, identifies, and controls reagent consumption and expiration dates



The Instrument Designed for You

FLEXIBLE PROGRAMMING

- Provides intuitive programming routines in multiple languages
- Each slide can be programmed to run a different protocol
- Provides reagent and slide loading maps
- A convenient run time display
- Programmable delayed run start

Maintain Profitability and High Quality

COST EFFECTIVE

- Cost efficient, minimal reagent usage
- Consumes less than 2 liters of buffer for a full 36 slide run
- Significant less waste output
- Separate containers for hazardous and non-hazardous waste to mitigate waste expense

Reproducible Brilliant Stains at Any Time

OPTIMIZED PROTOCOLS

- To further enhance throughput and the overall quality of the staining process, the DBS HighLighter[™] is delivered with optimized proprietary protocols for DBS antibodies and reagents
- Get access to unique protocols, customizable at any time

Fully Automated, High Quality, Efficient and Reliable



The Diagnostic BioSystems **HighLighter™**

SYSTEM ORDERING INFORMATION

DBS HIGHLIGHTER™ 360 AUTOMATED SLIDE STAINER

PART NUMBER: HL-AS360

The DBS HighLighter™ comes complete with the following components:

PART NUMBER	DESCRIPTION
HL-PC706B	2 x Reagent Rack, Modified
HL-PC326	Power Cord, Instrument
HL-324	USB Cable
HL-PC417A	2 Liter Bottle with Modified Cap and Fitting
HL-PC423	2 x 4 Liter Bottle with Modified Cap
HL-PC704	Label Printer
HL-PC716	Labels, Xylene Resistant
HL-PC717	Ribbon, Xylene Resistant
HL-PC331	Computer
HL-PC7331A	Computer Mouse
HL-PC720	Slide Tissue Area Map—2 Slides
HL-PC740	Alcohol Swabs—1 Set of Boxes

Experience the
DBS HighLighter™
in action!

For a demo, email us at:
HighLighter@dbiosys.com

HIGHLIGHTER™ DIMENSIONS:

36.5 in x 22 in x 24 in/93 cm x 56 cm x 61 cm (W x H x D)

COMPUTER CONFIGURATIONS:

Laptop Computer, Windows 10 Operating System,
Printer: Zebra TLP 3844-Z or Zebra GX430t,
4 USB 2.0 Communication Ports



United States (Headquarters)

6616 Owens Drive, Pleasanton, CA 94588, USA
Toll Free: (888) 896-3350, customersupport@dbiosys.com



Contact your distributor
for more information
and visit us at www.dbiosys.com

Europe

Germany
Lake Constance
Cell: +49 7557 929 3915
wolfgang.vogel@dbiosys.com

Mexico

Avenida Aztecas # 95, Casa 4
Colonia Pueblo de los Reyes
Alcaldia Coyoacan
CDMX, C.P. 04330
Cell: +525535709693
a.herrera@dbiosys.com

India

612, Eden Square
St. John's Road
Secunderabad – 500003
Cell: +91 9958293222
anandam@dbiosys.com

Canada

32 Lemsford Drive
Markham, ON L3S 4H4
Cell: +1 416 219 2035
hosna.mujadidi@dbiosys.com