


Immunohistochemistry Humidity Chamber

Diagnostic BioSystems Slide Master365TM is a compact and efficient platform that facilitates Humidity required for IHC staining. Appropriate maintenance of humidity along with flat platform is of extreme importance for consistent immunohistochemistry staining.

Capabilities

- ▶ Ensures flat platform.
- ▶ Maintains humidity.
- ▶ Ease to rinse and avoids majority of individual slides handling during staining.
- ▶ Two variants, 40(SM365-40) and 20(SM365-20) slide staining capacity.



Specifications	40 Slides	20 Slides
Color	Black and White	Black and White
Weight	2.5 kg	1.5 kg
Dimensions	470 x 370 x 50 mm	370 x 270 x 50 mm
Slide Capacity	40 to 50 Slides	20 to 30 Slides
Optimal Run Temperature	Room Temperature	Room Temperature

Headquarters

6616 Owens Drive
Pleasanton, CA 94588, USA
Toll Free: (888) 896-3350
customersupport@dbiosys.com



www.dbiosys.com