


VWR[®] Bottle Top Dispenser

- 01. ONE DISPENSER FOR ALL APPLICATIONS
- 02. FULLY AUTOCLAVABLE
- 03. DESIGNED WITH ERGONOMICS AND INTUITIVE HANDLING





VWR® BOTTLE TOP DISPENSER

BOTTLE TOP DISPENSERS WITH HIGH CHEMICAL COMPATIBILITY AND ISO 8655 STANDARD ACCURACY

- Designed with ergonomics and intuitive handling in mind, VWR® Bottle Top Dispensers exhibit excellent chemical compatibility and offer reliable precision in practical laboratory environments
- Fully autoclavable at 121°C and 15psi for 10 to 15 minutes
- Can be used with solvents (refer to chemical compatibility chart)

Description	Accuracy	Increments	Precision	Cat. No.
VWR® Bottle Top Dispenser 0.25-2.5 ml	± 0.6%	0.05 ml	±0.2%	76319-582
VWR® Bottle Top Dispenser 0.5-5 ml	± 0.6%	0.1 ml	±0.2%	76319-584
VWR® Bottle Top Dispenser 1-10 ml	± 0.6%	0.2 ml	±0.2%	76319-586
VWR® Bottle Top Dispenser 2.5-30 ml	± 0.6%	0.5 ml	±0.2%	76319-588
VWR® Bottle Top Dispenser 5-60 ml	± 0.6%	1.0 ml	±0.2%	76319-590
VWR® Bottle Top Dispenser 10-100 ml	± 0.6%	2.0 ml	±0.2%	76319-592



800 932 5000 | VWR.COM

Prices, product, and/or services details are current when published and subject to change without notice. | Certain products or services may be limited by federal, state, provincial, or local regulations. | VWR, part of Avantor, makes no claims or warranties concerning sustainable/green products. Any claims concerning sustainable/green products are the sole claims of the manufacturer and not those of VWR International, LLC and/or Avantor, Inc. or affiliates. All prices are in US dollars unless otherwise noted. Offers valid in US and Canada, void where prohibited by law or company policy, while supplies last. | Trademarks are owned by Avantor, Inc. or its affiliates, unless otherwise noted. | Visit vwr.com to view our privacy policy, trademark owners, and additional disclaimers. © 2020 Avantor, Inc. All rights reserved.